

KSI Sample Stages



To meet the demand of optimal throughput and to take full advantage of the capabilities of the KSI Scanning Acoustic Microscopes an optimized sample stage is often a necessity.

KSI offers a wide range of sample stages to fulfill these requirements.


Key Advantages:

- Accurate, fast and efficient loading and unloading of specimens
- Minimized setup times
- Optimal sample fixation



Sample stages are available for:

- JEDEC TRAY's
- Wafer up to 12"
- Solar Cells
- IGBT
- Universal sample stage
- Trough scan sample stage
- Special vacuum sample stage for smallest samples
- Prepared for multihead measurement



Need a special sample stage?

Other customized sample stages available on request.

Contact us!



KSI Germany

Scanning Acoustic Microscopy

a member of
IP-holding GmbH
Walther-Rathenau-Str.18
35745 Herborn
Germany

Phone: +49 2772 - 47 66 055
info@ksisam.com
www.ksisam.com
www.acoustic-microscopy.com