


From Eye to Insight

Leica
MICROSYSTEMS



Reproducible results for Cryo-TEM

EM GP2 - AUTOMATIC PLUNGE FREEZER FOR EM GRIDS

- > Sensor controlled single-sided, flat blotting
- > Temperature and humidity controlled environmental chamber
- > Control of secondary cryogen temperature
- > Cryogen gas liquefier for secondary cryogen filling
- > Dry nitrogen gas flow environment to protect sample after freezing
- > Safety alarms for humidity and secondary cryogen temperature control
- > Window heater, stereo microscope and LED illumination for full process visibility
- > Docking station for TEM cryo transfer holders

CONNECT
WITH US!



Leica Microsystems CMS GmbH | Ernst-Leitz-Straße 17–37 | 35578 Wetzlar
Tel. +49 (0) 6441 29-0 | F +49 (0) 6441 29-2599

<https://www2.leica-microsystems.com/em-gp2>